



Avenue de la Gare BP 5
34240 LAMALOU LES BAINS

Tél. 33 (0) 467 955 927
Fax 33 (0) 467 957 554

serviceclient@medwin.fr
www.medwin.fr



Easy Moov
Enteral Feeding
Pumps
Intuitively
fingertips...



- > Intuitive
- > Functional
- > Secure





A full range of Enteral Feeding Pump and accessories with the line JUMBO for hospital used and EASY MOOV for Ambulatory and homecare used.

The **Easy Moov** Enteral feeding pumps had been developed to answer to Ambulatory and Homecare treatments working in any position and situation.

Easy Moov pumps line is the most EASY enteral feeding pump concept turn to the patients, the really "INTUITIVE" one.

- Weight less than 500g including battery and charge device.
- Battery Ni MH with 24 hours life and automatic charge when connected.
- Noiseless and permanent running (no step by step rotation)
- Flow rate from 1ml to 400ml/h and dose from 1ml to 9999ml.
- All requested visual and audio alarms such as Free flow- low battery – air- no set- occlusion high and down - end of dose - empty bag...) with possibility to activate or not each alarm.
 - Feeding and alarm history covering 290 events with capability to print it.
 - Computer system to do quick pump check up and control.
 - Screw adaptor for I.V.pole and rucksack for adult and children.
 - Accuracy +/- 10% of selected flow rate.
 - Clear touch panel control with selected Multi-language.



Easy Moov able to be used in all patient situations.



Références	Description
PM4	Enteral feeding pump Easy Moov with dose and volume feature
NUD145M	Pump set for Easy Moov with screw connector for soft bag
NUD145MP	Pump set for Easy Moov /pediatric use
NUUD145M	Pump set for Easy Moov with double line and screw connectors for soft bags
NUD145MCE	Pump set for bottle with basket
NP1GBD145M	Pump set for Easy Moov with 500ml bag attached
NP2GBD145M	Pump set for Easy Moov with 1200ml bag attached
ENTP1250	Rucksack children for Easy Moov
ENTP1251	Rucksack adult for Easy Moov
ENTP0035	Screw adaptor for I.V pole for Easy Moov
ENTP0020	Screw adaptor for rail for Easy Moov
4908183	Electric Sector adaptor.